AQUACEL™
Hydrofibre™ Dressing with Strengthening Fibre / Hydrofibre™ Verband mit Verstärkenden Fasern / Medicazione in Hydrofibre™ con Fibra Rinforzate / Pansement Hydrofibre™ avec Fibres Renforcées / Apósito de Hydrofibre™ con fibra reforzante / Pensa Hydrofibre™ com Fibra de Reforço / Hydrofibre™ Verband mit Verstärkenden Vezeln / Hydrofibre™ Förband med Förstärkande Fibra / Bandage med forstärkede fibre / Hydrofibre™ Bandasje med Forsterkende Fibre / Hydrofibre™ Umbüddi með Styrkjandi þrárum / Εμπίδοχο Hydrofibre™ με ενισχυτικές ίνες / Hydrofibre™verband mit verstärkenden Fasern
INSTRUCTIONS FOR USE / GEBRAUCHSANLEITUNG / ISTRUZIONI PER L'USO / CONSEILS D'UTILISATION / INSTRUCCIONES DE USO / INSTRUÇÕES DE USO / GEBRUIKSINSTRUCTIES / BRUKSANVISNING / KÄTTÖOHJEET / BRUGSANVISNING / BRUKSANVISNING / LEIBBEINIGAR UM NOTKUN / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / تعليمات الإستعمال
CONvaTec
ENGLISH
PRODUCT DESCRIPTION AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre is a soft, sterile, non-woven ribbon dressing composed of sodium carboxy-methylcellulose and regenerated cellulose fibres for strengthening. This conformable and highly absorbent dressing absorbs wound fluid and transforms into a soft gel, which maintains a moist environment to support the body's healing process and aid in the removal of nonviable tissue from the wound (autolytic debridement), without damaging new formed tissue.
INDICATIONS Under medical supervision AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre may be used for the management of: <ul style="list-style-type: none">leg ulcers, pressure ulcers (Stage II-IV) and diabetic ulcers, surgical wounds (eg., post-operative, wounds left to heal by secondary intention and donor sites), partial thickness burns, traumatic wounds (eg., abrasions and lacerations), exudate absorption in oncology wounds (eg., fungating cutaneous tumours, cutaneous metastases and Kaposi's sarcomas).
CONTRAINDICATIONS AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.
PRECAUTIONS AND OBSERVATIONS <ul style="list-style-type: none">•Sterility is guaranteed unless the pouch is damaged or opened prior to use. This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use. •Appropriate supportive measures should be taken where indicated (eg., use of graduated compression bandaging in the management of venous leg ulcers or pressure relief measures in the management of pressure ulcers/sores). •The control of blood glucose, as well as appropriate supportive measures, should be provided with diabetic foot ulcers. •Infection is not a contraindication to the use of AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre. Should infection develop during the use of the dressing, antibiotic therapy should be initiated, as clinically indicated, by a health care professional. •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre can facilitate the control of minor bleeding. •This ribbon product can be cut, if required prior to use, across the narrower dimension of the ribbon. Cutting along the longer dimension of the ribbon should not be performed. •If removing the dressing is difficult, the dressing should be fully saturated with sterile saline or water and removed slowly. •Because AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre provides a moist environment that supports the growth of new blood vessels, occasionally the delicate newly formed blood vessels may produce a blood stained wound fluid. •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre is not intended for use as a surgical sponge. •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre is not intended for use within internal body cavities or within closed wounds.
DIRECTIONS FOR USE For Exuding Wounds <ul style="list-style-type: none">•Before applying the dressing, cleanse the wound area with an appropriate wound cleanser. •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre should overlap at least 1 cm onto the skin surrounding the wound. •When using AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre in deep cavity wounds, leave at least 2.5 cm outside the wound for easy removal. Only pack directly into the wound to a depth of 1 cm. AQUACEL™ Dressing with Strengthening Fibre will expand to fill the wound spaces on contact with wound fluid. •Apply the dressing to the wound and cover with a moisture retentive dressing (eg., DuoDERM™ Extra Thin), or other appropriate secondary dressing. See relevant individual product package inserts for complete instructions for use. •All wounds should be inspected frequently. Remove the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre when clinically indicated (i.e., leakage, excessive bleeding, or suspicion of infection). •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre can be left in place for up to 7 days, where clinically indicated, (except in burns where it can be left in place for up to 14 days).
FOR DRY WOUNDS In addition to the directions for use set forth above: <ul style="list-style-type: none">•Place the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre on the wound and wet with sterile water or saline over the wound area only. •Cover the dressing with a moisture retentive dressing such as DuoDERM™ Extra Thin to avoid drying out of the dressing and subsequent dressing adherence to the wound.
FOR PARTIAL THICKNESS BURNS •Before applying the dressing, cleanse the wound area with an appropriate wound cleanser. <ul style="list-style-type: none">•The AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre should overlap at least 5 cm onto the skin surrounding the burn or other adjacent AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre, as the dressing will shrink as it absorbs the exudate. •The AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre should be covered with a sterile loosely woven gauze pad and appropriately secured. In the immediate post burn period (up to 4 days) large volumes of wound exudate may require the use of a saturated dressing be removed and replaced with new AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre. •Remove the cover dressing periodically and inspect the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre while it remains in place on the burn. •If complete adherence to the wound bed of the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre is a desired characteristic. Non adherence of areas of the dressing may indicate deepening of the wound or infection. Areas of the dressing may be cut away to facilitate assessment. The exposed areas of the dressng may be treated appropriately. •As the burn wound re-epithelializes, the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre will detach or be easily removed. •For partial thickness burns, AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre may be left in place for up to 14 days or until clinically indicated. •If complete adherence to the wound bed of the dressing may be necessary.
Discard any unused portion of the dressing
If the immediate product packaging is damaged, do not use. Store in a cool dry place.
If further information or guidance is needed, please contact ConvaTec Professional Services.

CONvaTec Limited First Avenue, Deeside Industrial Park, Deeside, Flintshire, CH5 2NU, UK
CONvaTec Inc. Greensboro, NC 27409 USA 1-800-422-8811
www.convatec.com
© 2018 ConvaTec Inc. */™ indicates a trademark of ConvaTec Inc.

GERMAN
PRODUCT DESCRIPTION AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre is a soft, sterile, non-woven ribbon dressing composed of sodium carboxy-methylcellulose and regenerated cellulose fibres for strengthening. This conformable and highly absorbent dressing absorbs wound fluid and transforms into a soft gel, which maintains a moist environment to support the body's healing process and aid in the removal of nonviable tissue from the wound (autolytic debridement), without damaging newly formed tissue.
INDICATIONS Under medical supervision AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre may be used for the management of: <ul style="list-style-type: none">leg ulcers, pressure ulcers (Stage II-IV) and diabetic ulcers, surgical wounds (eg., post-operative, wounds left to heal by secondary intention and donor sites), partial thickness burns, traumatic wounds (eg., abrasions and lacerations), exudate absorption in oncology wounds (eg., fungating cutaneous tumours, cutaneous metastases and Kaposi's sarcomas).
CONTRAINDICATIONS AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or its components.
PRECAUTIONS AND OBSERVATIONS <ul style="list-style-type: none">•Sterility is guaranteed unless the pouch is damaged or opened prior to use. This device is for single-use only and should not be re-used. Re-use may lead to increased risk of infection or cross contamination. Physical properties of the device may no longer be optimal for intended use. •Appropriate supportive measures should be taken where indicated (eg., use of graduated compression bandaging in the management of venous leg ulcers or pressure relief measures in the management of pressure ulcers/sores). •The control of blood glucose, as well as appropriate supportive measures, should be provided with diabetic foot ulcers. •Infection is not a contraindication to the use of AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre. Should infection develop during the use of the dressing, antibiotic therapy should be initiated, as clinically indicated, by a health care professional. •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre can facilitate the control of minor bleeding. •This ribbon product can be cut, if required prior to use, across the narrower dimension of the ribbon. Cutting along the longer dimension of the ribbon should not be performed. •If removing the dressing is difficult, the dressing should be fully saturated with sterile saline or water and removed slowly. •Because AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre provides a moist environment that supports the growth of new blood vessels, occasionally the delicate newly formed blood vessels may produce a blood stained wound fluid. •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre is not intended for use as a surgical sponge. •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre is not intended for use within internal body cavities or within closed wounds.
DIRECTIONS FOR USE For Exuding Wounds <ul style="list-style-type: none">•Before applying the dressing, cleanse the wound area with an appropriate wound cleanser. •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre should overlap at least 1 cm onto the skin surrounding the wound. •When using AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre in deep cavity wounds, leave at least 2.5 cm outside the wound for easy removal. Only pack directly into the wound to a depth of 1 cm. AQUACEL™ Dressing with Strengthening Fibre will expand to fill the wound spaces on contact with wound fluid. •Apply the dressing to the wound and cover with a moisture retentive dressing (eg., DuoDERM™ Extra Thin), or other appropriate secondary dressing. See relevant individual product package inserts for complete instructions for use. •All wounds should be inspected frequently. Remove the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre when clinically indicated (i.e., leakage, excessive bleeding, or suspicion of infection). •AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre can be left in place for up to 7 days, where clinically indicated, (except in burns where it can be left in place for up to 14 days).
FOR DRY WOUNDS In addition to the directions for use set forth above: <ul style="list-style-type: none">•Place the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre on the wound and wet with sterile water or saline over the wound area only. •Cover the dressing with a moisture retentive dressing such as DuoDERM™ Extra Thin to avoid drying out of the dressing and subsequent dressing adherence to the wound.
FOR PARTIAL THICKNESS BURNS •Before applying the dressing, cleanse the wound area with an appropriate wound cleanser. <ul style="list-style-type: none">•The AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre should overlap at least 5 cm onto the skin surrounding the burn or other adjacent AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre, as the dressing will shrink as it absorbs the exudate. •The AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre should be covered with a sterile loosely woven gauze pad and appropriately secured. In the immediate post burn period (up to 4 days) large volumes of wound exudate may require the use of a saturated dressing be removed and replaced with new AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre. •Remove the cover dressing periodically and inspect the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre while it remains in place on the burn. •If complete adherence to the wound bed of the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre is a desired characteristic. Non adherence of areas of the dressing may indicate deepening of the wound or infection. Areas of the dressing may be cut away to facilitate assessment. The exposed areas of the dressing may be treated appropriately. •As the burn wound re-epithelializes, the AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre will detach or be easily removed. •For partial thickness burns, AQUACEL™ Hydrofibre™ Dressing with Strengthening Fibre may be left in place for up to 14 days or until clinically indicated. •If complete adherence to the wound bed of the dressing may be necessary.
Discard any unused portion of the dressing
If the immediate product packaging is damaged, do not use. Store in a cool dry place.
If further information or guidance is needed, please contact ConvaTec Professional Services.

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© 2018 ConvaTec Inc. */™ kenzeichnen Marken von ConvaTec Inc.