# Tutoplast® Spongiosa

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Tutoplast® Spongiosa

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## 4.3 Contraindications

None known

### 4.4 Special warnings and precautions for use

- Strict assessment of the indications must be performed in the following cases prior to use of the graft:
  - For implantation in necrotic host sites
  - For implantation in regions with poor blood supply
  - For implantation in host sites with an active or latent infection
  - In the case of disorders or conditions which could have a negative impact on the healing rate
- Based on their medical history, as well as further tests, all tissue donors are examined for possible grounds for exclusion. Before the donor tissue can be approved for processing, serological tests are carried out for possible infectious agents. This analysis includes tests for HTLV-I/-II, HIV-1/-2, hepatitis B/C and Treponema pallidum, among others.
- The Tutoplast® process is effective against all known pathogenic species. However, as with all biological products, the transmission of infectious diseases can be ruled out only to the extent of current knowledge.
- The graft remains sterile as long as the packaging is undamaged. If the sterile packaging becomes damaged, the product must be disposed of appropriately. Do not resterilise the product! If the graft is contaminated in the course of a surgical intervention, it must be disposed of appropriately.
- When using the graft, note that the mechanical loads vary depending on the implant site and the physical stability of the product must be sufficient for the respective conditions.
- When customising the graft for the individual patient, do not allow it to be damaged by the heat generated. The structural integrity must be preserved.
- The graft must be implanted (ideally using the "press-fit" technique) and fixated in such a way that it remains in the intended site and its incorporation is facilitated.
- When inserting the blocks, it is important to avoid exerting loads exceeding the physical stability of the graft under all circumstances, so as to preserve the structural integrity of the graft. This is particularly important for applications in spinal surgery.
- In particular in the case of patient-specific graft blocks, which are manufactured according to the specifications of the attending physician, the type of fixation used on thin regions must be considered particularly carefully.
- Postoperative loads exceeding the physical stability of the graft must be accommodated in additional measures.
- Tutoplast® Cancellous is intended for single use only. Unused material must be disposed of appropriately.

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## 4.5 Interaction with other medicinal products and other forms of interaction

None known

### 4.6 Pregnancy and lactation

No clinical data concerning the use of the product during pregnancy and when breastfeeding are available. Based on the physical mechanism of action of the graft, undesirable effects during pregnancy or on the health of the newborn child are not to be expected and have not been determined.

### 4.7 Effects on ability to drive and use machines

Not relevant

### 4.8 Undesirable effects

The following classifications are used as the basis for the frequency of adverse reactions:

Very common (>1/10)

Common (>1/100 to <1/10)

Uncommon (>1/1.000 to <1/100)

Rare (>1/10.000 to <1/1.000)

Very rare (<1/10.000)

Not known (frequency cannot be estimated from the available data)

System organ class	Frequency	Adverse effect
Immune system disorders	Not known	Graft rejection
General disorders and administration site conditions	Not known	Reaction at the site of implantation <p>Dehiscence</p>
Injury, poisoning and procedural complications	Not known	Graft failure

As with every surgical procedure, there is the possibility of infection and wound complications with associated symptoms (swelling, redness, inflammation and pain) or other reaction due to the procedure itself.

Reporting of suspected adverse reactions

Reporting suspected adverse reactions of the medicinal product after its approval is important. It allows continued monitoring of the benefit/risk balance of the medicinal product.

- Healthcare professionals within Germany are asked to report any suspected adverse reactions to the Bundesinstitut für Impfstoffe und biomedizinische Arzneimittel (Federal Institute

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- Healthcare professionals outside Germany are asked to report any suspected adverse reactions to the manufacturer:

Tutogen Medical GmbH, Industriestrasse 6, 91077 Neunkirchen am Brand, Germany Tel.: +49 9134-9988-0, Fax: +49 9134-9988-99, vigilance@rtix.com

National legislation may stipulate additional requirements for reporting to local authorities.

Patients must be informed to contact their doctor or medical professional if they notice any side effects. This also applies to side effects that are not stated in this leaflet and summary of product characteristics. Patients can also report side effects directly to the Paul Ehrlich Institute. By reporting side effects, patients can help provide more information about the safety of this medicinal product.

## 4.9 Overdose

Not relevant

## 5. Pharmacological properties

### 5.1 Pharmacodynamic properties

Pharmacotheapeutic group: Allograft

*Mechanism of action of the medicinal product*

The medicinal product is primarily a physical act, acting as a guide rail for the regeneration of osseous tissue. As it can be absorbed by osteoclasts, it is integrated into the autologous remodelling process. The regeneration of autologous structures (ossification/osteoegenesis) occurs in parallel to the resorption of the allograft. Both processes are particularly dependent on the host site's ability to respond, its blood supply and the size of the graft itself. The implantation of the graft and its function as a guide rail facilitate the naturally occurring process of bone remodelling at the site of the defect. The ultimate goal is complete conversion of the graft into new autologous tissue.

### 5.2 Pharmacokinetic properties

Tutoplast® Cancellous is primarily a physical act, acting as a guide rail for the regeneration of osseous tissue. The process of cell colonisation, remodelling and resorption generally begins after 2 days. The conversion process depends on the size of the graft and the host site's ability to respond. As such, it can take a considerable period of time to complete (several months or even years).

## 5.3 Preclinical safety data

Studies on animals and *in vitro* studies have shown that the graft is biocompatible, is incorporated quickly and maintains the space until replaced by autologous tissue. As such, the preclinical data do not indicate any risks to humans and confirm that the product can be used safely.

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None known

**6.3 Shelf life**

Tutoplast® Cancellous has a shelf life of 5 years providing the original packaging remains intact. The graft must not be used after the expiry date shown on the packaging.

**6.4 Special precautions for storage**

The graft should be stored in a clean, dry place protected from exposure to direct sunlight and at no higher than 30°C. The product must not be frozen.

**6.5 Nature and contents of container**

The graft is kept sterile in a double-layered primary packaging. Blocks are additionally packaged in a tubular bag, granulate in a canister. The packaging contains either a pre-shaped block or the quantity of granulate stated on the external packaging.

Shape, size	Packaging unit
Half Femoral head	Half Femoral head
Customised Block	Min. 5 mm x 1 mm x 1 mm Max. 60 mm x 30 mm x 30 mm
Block	8 mm x 8 mm x 8 mm, 20 mm x 20 mm x 10 mm, 10 mm x 10 mm x 20 mm, 10 mm x 20 mm x 20 mm, 18 mm x 15 mm x 6 mm, 20 mm x 10 mm x 10 mm, 20 mm x 20 mm x 30 mm, 30 mm x 20 mm x 12 mm, 30 mm x 10 mm x 10 mm, 30 mm x 30 mm x 12 mm
Block-P	15 mm x 10 mm x 9 mm, 15 mm x 15 mm x 9 mm
Dowel	Ø 7 mm, length 14 mm - 18 mm, Ø 10 mm, length 16 mm - 20 mm, Ø 12 mm, length 16 mm - 20 mm, Ø 14 mm, length 16 mm - 20 mm, Ø 16 mm, length 16 mm - 20 mm

Wedge	35 mm x 22 mm x 5 mm, 35 mm x 22 mm x 8 mm, 35 mm x 22 mm x 10 mm, 35 mm x 25 mm x 10 mm, 35 mm x 25 mm x 12 mm	1 piece
Disc	Ø 3 mm, 5 mm, 8 mm, 10 mm, 15 mm	1 piece
Iliac Crest Wedge	7 mm x 25 mm, 10 mm x 25 mm, 12 mm x 25 mm, 13 mm x 30 mm, 16 mm x 25 mm, 20 mm x 25 mm, 26 mm x 25 mm, 40 mm x 25 mm	1 piece
Cortico-cancellous Chips	Ø 0,25 mm - 1 mm	0,5 cm³, 1 cm³, 2 cm³
	Ø 1 mm - 2 mm	0,5 cm³, 1 cm³, 2 cm³
Chips	Ø 2 mm - 4 mm	3 cm³, 5 cm³, 10 cm³, 15 cm³, 30 cm³
	Ø 4 mm - 10 mm	10 cm³, 15 cm³, 30 cm³
	Ø 0,25 mm - 1 mm	0,5 cm³, 1 cm³
	Ø 1 mm - 2 mm	0,5 cm³, 1 cm³, 2 cm³, 3 cm³, 5 cm³, 10 cm³, 15 cm³
Chips	Ø 1 mm - 4 mm	15 cm³, 30 cm³
	Ø 1 mm - 10 mm	5 cm³, 10 cm³, 15 cm³, 20 cm³, 30 cm³
	Ø 2 mm - 4 mm	3 cm³, 5 cm³, 10 cm³, 15 cm³, 20 cm³, 25 cm³, 30 cm³
	Ø 4 mm - 10 mm	5 cm³, 10 cm³, 15 cm³, 20 cm³, 25 cm³, 30 cm³

Note: It is possible that not all packaging sizes will be available on the market.

**6.6 Special precautions for disposal and other handling**

**Information on disposal of remnants**  
If the packaging is opened and the graft not used, the cancellous bone must not be resterilised.

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Remnants must be disposed of in accordance with the pertinent national regulations. The contents of the packaging are exclusively intended for use in **one** patient.

**Special information on handling**

Please observe the general principles of aseptic handling when using Tutoplast® Cancellous.

- The graft must be completely rehydrated in normal saline solution or Ringer's solution prior to its use.
- If the graft is to be customised before use, this should only be done once it is completely rehydrated.
- Any dust which forms must be thoroughly rinsed away with normal saline solution or Ringer's solution.
- Autologous blood, blood components or bone marrow aspirate may only be added after rehydration and immediately prior to the product's implantation.

**a) Standard rehydration for all forms**

- Immerse the graft completely in a normal saline solution (or Ringer's solution) in a sterile container – the deeper the water, the better.
- Weight blocks down with a sterile instrument to prevent them from floating to the surface.
- The rehydration time should be at least 30 minutes.
- During the rehydration process, agitate the cancellous bone and/or shake the container repeatedly.
- Rehydration is complete when no more air bubbles escape the cancellous bone when it is agitated and the bone remains on the bottom of the container without being weighed down instead of floating to the surface.
- Store the graft in the rehydration solution until ready for implantation.

**b) Vacuum rehydration for all forms** (vacuum rehydration permits quicker, more efficient rehydration and is thus particularly recommended for blocks)

- The graft is placed in an appropriately sized, sterile, disposable syringe.
- Draw up normal saline solution (or Ringer's solution) with the filled syringe until the graft is completely covered.
- Turn the syringe vertical and eject all the air.
- Seal off the opening of the syringe so that it is airtight (using a gloved finger, for example).
- Pull the plunger down with force. Holding the plunger retracted, shake the syringe back and forth repeatedly. Air bubbles can be seen escaping the cancellous bone. Repeat the procedure a number of times.
- When correctly rehydrated, the bone will settle at the bottom of the syringe. When this occurs, repeat the procedure at least once more.
- Store the graft in the rehydration solution until ready for implantation.

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07/2024 "11"

**11. Sales restrictions (DE)**

Only available on prescription  
Notes:  
Each graft is supplied with stickers for documentation in the patient's file. Tutoplast® is a registered trademark of Tutogen Medical GmbH.

**Tutogen Medical GmbH**  
Industriestraße 6  
91077 Neunkirchen am Brand  
Germany  
Tel.: +49 (0)9134 9988-0  
Fax: +49 (0)9134 9988-99  
tutogen@trix.com

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