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Instructions to User

This Fingertip Pulse Oximeter is a kind of innovated medical device with non-invasive and continuous features for artery SpO₂ and PR detection. Being portable, it is able to measure SpO₂ and PR values quickly and precisely. This product is medical device, which can be used repeatedly.

General Description


Haemoglobin Saturation is the percentage between the capacity of Oxyhemoglobin (HbO₂) that compounded with oxygen and that of all combinative haemoglobin (Hb) in blood. In other words it is the saturation of Oxyhemoglobin in blood. It is a very important physiological parameter for Respiratory and Circulation Systems. Many respiratory diseases could reduce haemoglobin saturation in human blood. Moreover, factors such as Automatic Organic Regulation Malfunction caused by anaesthesia, trauma resulted from major operation and some medical examination can also cause problems in oxygen supply, which might reduce human haemoglobin saturation. As a result, such symptoms as megrim, vomiting and asthenia might appear to patients. Hence, it is very important to know haemoglobin saturation of patient timely in clinical medical aspects. The Fingertip pulse oximeter features in small volume, low power consumption, convenient operation and portability. It is only necessary for patient to put one of finger into fingertip photoelectric sensor for diagnosis, and the display screen will directly show the measured value of hemoglobin saturation. It has been proved in clinical experiments that it possesses rather high precision and repeatability.

Measurement principle

The principle of the oximeter is as follows: An experience formula of data process is established by exerting Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (R Hb) and Oxyhemoglobin (O₂ Hb) in glow and near-infrared zones. Operation principle of the instrument is to combine Photoelectric Oxyhemoglobin Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that two lights with different wavelength (660nm glow and 905nm near-infrared light) can be focused onto human nail through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element, information acquired through which will be shown on two groups of LEDs through process in electronic circuits and microprocessor.

Diagram of Operating Principle

1. Infrared-ray receiving tube
2. Infrared-ray transmitting tube

Precautions for Use

1. Do not use the pulse oximeter together with MRI or CT equipment.
2. Explosion hazard: Do not use the pulse oximeter in an explosive atmosphere.
3. The pulse oximeter is intended only as an adjunct in patient assessment. Doctors should make the diagnosis in conjunction with clinical manifestation and symptoms.
4. Check the pulse oximeter sensor application site frequently to make sure that the circulation and skin integrity of patient are under good condition.
5. Do not stretch the adhesive tape while applying the pulse oximeter sensor. This may cause inaccurate reading or skin blisters.
6. Please read the manual carefully before operation.
7. Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status and correct alignment at least every 2 hours.
8. Inaccurate measurements may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid.
9. Significant levels of dysfunctional hemoglobins (such as carboxyl hemoglobin or methemoglobin) may cause inaccurate reading.
10. Intravascular dyes such as indocyanine green or methylene blue may cause inaccurate reading.
11. SpO₂ measurements may be adversely affected in the presence of high ambient light. Please shield the sensor area (e.g. with your hand or a cloth) if it is necessary.

12. Unexpected action may cause inaccurate reading.
 13. Medical signal with high frequency or interference caused by a defibrillator may lead to inaccurate reading.
 14. Venous pulsations may cause inaccurate reading.
 15. It may cause inaccurate reading when the positions of sensor and blood pressure cuff are on the same arterial catheter or intravascular line.
 16. Hypotension, severe vasoconstriction, severe anemia, or hypothermia may cause inaccurate readings.
 17. It may cause inaccurate reading by giving use of cardiogenic to patient after his cardiac arrest or when he is in quiver.
 18. Bright nail or painted nail may cause inaccurate SpO2 reading.
- Follow local ordinances and recycling instructions regarding to disposal or recycling of the device and device components.

Scope of Application/Intended Use

The fingertip pulse oximeter can be used to measure the human haemoglobin saturation and heart rate through finger. It can be used in hospitals, families, schools and medical centers.

Contraindications: not found

NOTE:

1. The image in the instruction may have slight differences with the actual instruments.
2. Technical parameters and appearance change without prior notice. Product include: main machine and SpO2 sensor.

Features

- OLED Display: Product adopts double color OLED display, can show the six different display mode.
- Low power consumption, continuously work for more than six hours with one rechargeable lithium battery.
- Battery power indicator.
- In the absence of signals, the product will power-off automatically after 8 seconds.
- Small in volume, light in weight, and convenient to carry.

Operation instructions

1. Plug one finger into rubber hole of the Oximeter (it is best to plug the finger thoroughly) before releasing the clamp with the nail upwards.
2. Press button on the front panel.
3. Don't tremble your finger when the Oximeter is working. Your body is not recommended on moving status.
4. Press the button on the front panel, if you want to change the display direction.
5. Read the relevant data from the display screen.
6. Please charge the battery in time when OLED indicates that the battery are in low power.

When plugging your finger into the Oximeter, your nail surface must be upward.

Declaration: Please use the medical alcohol to clean the rubber before each test and clean the tested finger with alcohol before and after the test. (The rubber inside of the Oximeter adopts medical rubber, which has no toxin, no harm, and brings no side effect such as allergy to the skin.)

Brief description of front panel

Characteristics

1. OLED Display

2. Charging jack

3. Start / mode Key

OLED Display Modes

Keys function description: In standby mode, start the device into the working state. Push down the button under working state to change the display mode. The device picture only for your reference use, specific in kind prevail.

Parameter setup

Press start button (>3s), into parameter setup as menu 1.

1. In menu 1, when the "*" signal is shown on the "Alm Setup", press the button (>3s) and enter into the menu 2. Press the button (<1s) can select item, then press button (>3s) to set the on/off for prompt, beep, demo and screen brightness adjustment (optional 1, 2, 3 and 4). When the "*" signal is shown on the "Restore", press the button (>3s) and all settings are back to factory settings.

Attention

Using devices outdoors or under strong light, please adjust the screen brightness to a higher proper level for observation. It's better for user to choose a lower brightness to conserve battery power.

2. In menu 2, when the "*" signal is shown on the "Sounds Setup", press the button (>3s) and enter into the menu 1. Press the button (<1s) can select item, then press button (>3s) to setup data. Choose + or - to plus or minus values.

Settings		Settings	
Alm Setup	*	Sounds Setup	*
Alm	off	SpO2 Alm Hi	100
Beep	off	SpO2 Alm Lo	100
Demo	off	PR Alm Hi	100
Restore	ok	PR Alm Lo	100
Brightness	3	+/-	+
Exit		Exit	
Menu 1		Menu 2	

Detailed Description of Product Functions

1. Display type: OLED-Display
2. SpO2 Measurement range: 70%—99%
Accuracy: $\pm 2\%$ on the SpO2 80%-99%:
 $\pm 3\%$ on the SpO2 70%-79%
Below 70% no requirement
Resolution: $\pm 1\%$
3. PR Measurement range: 30BPM-240BPM
Accuracy $\pm 1\text{BPM}$ or $\pm 1\%$ (the larger one)
4. Parameters of LED probe
Red light: Wavelength $660 \pm 2\text{nm}$, Radiation power 1.8mW
IR light: Wavelength $905 \pm 10\text{nm}$, Radiation power 2.0mW
5. Power: Rechargeable lithium battery
6. Automatic standby: The product shuts off by itself when no finger is inserted for about 8 seconds.
7. Dimension: 44mm x 28,3mm x 26,5mm
8. Operation environment: Operation temperature: $5^{\circ}\text{C} - 40^{\circ}\text{C}$
Storage temperature: $-10^{\circ}\text{C} - 40^{\circ}\text{C}$
Ambient humidity: 15%-80% in operation
10%-80% in storage
Air pressure: 86kPa-106kPa
9. Declaration: EMC of this product comply with IEC60601-1-2 standard.
10. Measurement Performance in Low Perfusion Condition: Required test equipment (BIO-TEK INDEX

Pulse Oximeter tester) could measure the available pulse wave with the amplitude of 6% of the simulation pulse wave amplitude.

11. Interference Resistance Capacity against ambient light: Device works normally when BIO-TEK INDEX Oximeter tester exerts interfering signal test.

Classification

1. Management Class for Medical Devices: II Equipment
2. Anti-electric Shock Type: Internally powered equipment
3. Anti-electric Shock Degree: Type BF equipment

Maintenance and Preservation

1. Charge batteries timely when the battery power is too low.
2. Clean the surface of fingertip oximeter before it is used to diagnose patients.
3. It is recommended to preserve the product in -10 - 40°C (14 - 104°F) and humidity is 10%-80%.
4. It is recommended that the product should be kept dry any time. A wet ambience might affect its lifetime and even damage the product.

Product Accessories

1. User manual
2. Charging cable

Guidance and Manufacturer's Declaration - Electromagnetic Radiation for Other Equipments and Systems

The Pulse Oximeter is designed to be used in specified electromagnetic environment. Users of the Pulse Oximeter must use it in the following environments.		
Radiation Test	Compliance	Electromagnetic environment guidance
RF-interference CISPR 11	Group 1	RF signal of Pulse Oximeter is simply created by its internal function. Therefore, its RF interference is very low and is not likely to cause any interference to nearby electronic equipment.
RF-interference CISPR 11	Class B	The Pulse Oximeter applies to all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Troubleshooting

Problem	Possible Reason	Solution
The SpO2 and Pulse Rate can not be shown normally	1. Finger is not properly positioned 2. Patient's SpO2 is too low to be measured.	1. Place the finger properly and try again. 2. Try again; If you are sure the device works all right, go to a hospital timely for exact diagnosis.
The SpO2 and Pulse Rate are shown unsteadily	1. Finger is not placed inside deep enough. 2. Finger is trembling or the patient is moving.	1. Place the finger properly and try again. 2. Please remain at rest.
The device can not be turned on	1. Inadequate power or power off. 2. Oximeter might be damaged.	1. Charge battery timely. 2. Please contact the local customer service.
Display is suddenly off	1. The device automatically powers off when it gets no signal within 8 seconds. 2. Inadequate power.	1. Normal. 2. Charge battery timely.

Symbols and Definitions

Symbol	Description	Symbol	Description
	Type BF applied part		European Representative
	Swiss representation		
	WEEE (2002/96/EC)		Serial number
	Refer to instruction manual		Lot number

	Caution		Date of manufacture
IP22	International Protection degree		Manufacturee
	Humidity range		Temperature range
	Keep dry		Avoid sunlight
	Product CE certification		ON/OFF

Reserves the right to technical change appearance, our products are subject to change without prior notice.

Novidion GmbH
 Fuggerstr. 30 • 51149 Colonia • Alemania
 Tel.: +49 (0) 2203 / 9885 200
 Fax: +49 (0) 2203 / 9885 206
 www.pulox.de • Mail: info@novidion.de

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Xuzhou Yongkang Electronic Science Technology Co., Ltd
Prolinx GmbH
MedNet SWISS GmbH
Novidion GmbH, Fuggerstr. 30, 51149 Köln