INHALATOR IN-100

Bedienungsanleitung
Instruction Manual
Istruzioni per l'uso
Mode d'emploi
Instrucciones de uso
Gebruiksaanwijzing
Instrukcja obsługi
Bruksanvisningar

www.pulox.de

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Dear Customer,

we are pleased that you have chosen a product from our range. Please read these operating instructions carefully and keep them accessible for other users. Please observe the notes. We do not accept any responsibility for damage resulting from failure to observe the instructions for use. We reserve the right to update and modify the instructions for use without prior notice.

With kind regards

Your Novidion Team

Notice:

Please place the masks in warm water for a few minutes before first use. This will give them an optimal fit.

1. Safety

Please read the operating instructions carefully to ensure safe use.

Warning: Indicates that a procedure is dangerous or unsafe and that death, personal injury or property damage may occur if the operation is continued.

Attention: Indicates important notes, instructions or explanations for optimal use.

Warning

- Comply with the type of medication, dosage and method of use prescribed by your doctor. Otherwise, there may be a worsening

of symptoms.

- Follow the operating modes specified in the instructions for use, as otherwise malfunctions may occur.
- The device accessories are designed for single use; do not use the accessories repeatedly, as this may lead to cross-infections.
- When using the device for the first time or if the medication cup has not been used for a long time, clean and disinfect cup and mask. Otherwise, bacterial infections may occur.
- Each user must use their own accessories, otherwise cross-infection may occur.
- Clean the accessories after disinfection, otherwise the patient could inhale the remaining disinfectant, which could lead to a worsening of symptoms.
- Do not reuse used medication; use new medication for each application. Otherwise, the patient could become infected with bacteria, which could lead to a worsening of symptoms.
- Do not use the device to inhale water, otherwise there may be a worsening of symptoms.
- Do not use the device at ambient temperatures above 40°C. There may be injury to the nasal mucosa or failure of the device.
- Do not rinse the main unit with water, immerse it in water, or store it in a humid environment. It may cause device failure.
- Clean the device after use and dry it immediately after cleaning. Otherwise, the patient may be infected by bacteria.
- Keep the device out of the reach of children, as there is a risk of small parts being swallowed.
- Do not use the device in the presence of combustible or explosive gases or anesthetic mixtures. Otherwise, injury may result.
- Avoid children placing dangerous parts such as cables around their necks, as this can cause suffocation.

Attention:

- If the device does not automatically switch off when there is a lack

of liquid, please press the “ON/OFF” button immediately to avoid damaging the device. See chapter 6 Troubleshooting.

- Clean the medication cup after each use. Otherwise, the device will not work correctly.
- The recommended maximum liquid capacity is 10ml, and the minimum capacity is 2ml.
- When the device is used at its maximum capacity and the temperature of the medication in the cup is higher than the ambient temperature during normal use, the maximum temperature should be below 40°C.
- When cleaning the medication container, do not place the device directly under a water jet, as water may enter the device.
- Do not use the device near high-frequency electromagnetic transmitters or other high-frequency electronic devices.
- Hold the device as vertically as possible during use.
- Do not drop the main unit or medication container and do not subject them to strong impacts.
- Do not touch the metal mesh of the nebulizer plate with a cotton swab or sharp objects. Otherwise, it may cause functional damage to the device.
- Patients who have sensitive areas with bruises, burns, inflammation, and facial/oral injuries should avoid use. If any discomfort occurs during use, please stop using it immediately and consult a doctor.
- Commonly used medications include mucolytics, bronchodilators, and antibiotics, such as terbutaline sulfate solution for nebulization and ipratropium bromide solution for inhalation. The active ingredient is water-soluble, non-irritating, non-toxic and does not trigger any allergic reactions. The pH value is almost neutral, it can adapt to the colloid osmotic pressure of the tissue and has a good nebulization effect and stability.
- The face mask is made of softener-free PVC material. Clinical medical personnel should follow the instructions of the

medication to determine whether the PVC accessories can be used.

- The product can be reused. Please disinfect it immediately after use, otherwise a bacterial infection may occur.
- This product should be purchased and used under a doctor's guidance.
- Do not use highly concentrated medical solutions. Please consult your doctor for specific types of nebulization medications and follow the doctor's recommendations to the user.
- Charge the device when the battery power seems to be low.
- If the device is not used for a long period of time, please recharge it regularly.
- Ensure that a legal guardian is present when children use it.
- Do not store or transport the device with medication in the medication container.
- Disposal of old devices and accessories must be carried out in accordance with local authority regulations.
- The use of this product differs from moistening the larynx and nasal mucosa.
- This product cannot be used in respiratory anesthesia and ventilation systems.
- The product life is 3 years (excluding consumables).
- The accessories supplied with the device are disposable and processed by ethylene oxide sterilization. Please check the packaging before use. Do not use the accessories if the packaging is damaged and contact the supplier.
- We can provide schematics, component lists, and other information needed for servicing. To do so, please contact the supplier.
- Date of manufacture: see label

2. Overview

2.1 Function and use

Intended use:

The product is intended for the nebulization of liquid medications for inhalation by the patient. The device can be used on children and adults at home, in the hospital and in outpatient aftercare.

Indication: Nebulizer therapy

Patient group: Adults and children

Intended users: Healthcare professionals and patients under their direction

Contraindications: The product must not be used together with pentamidine.

2.2 Features

Power supply: DC 5 V or 3.7 V rechargeable lithium battery

Input power: <5 VA

Atomization rate in level I: ≥ 0.25 ml/min

Atomization rate in level II: ≥ 0.15 ml/min

Noise level: ≤ 50 dB

Equivalent volume-particle size distribution: The volume distribution ratio of small atomized particles (diameter $< 5 \mu\text{m}$) is at least 60%

Type of protection against electric shock: Class II device, internally powered device

Degree of protection against electric shock: Type BF applied part

Degree of protection against ingress of liquid: IP22

Note: Please purchase a medical-grade power adapter from a qualified manufacturer (input: AC 100-240V, 50-60Hz; output: DC 5V, 1.0A).

2.3 Operating environment

Temperature: $5^{\circ}\text{C} \sim 40^{\circ}\text{C}$

Air humidity: 15 % $\sim$ 90 %.

Air pressure: 700 hPa $\sim$ 1060 hPa

Caution: This product is not suitable for use in environments with strong electromagnetic interference (e.g. various medium/high

frequency therapeutic instruments, transformers, large electrical cabinets, radio and television transmission towers, other high frequency transmission devices and other electrical equipment or medical devices that may generate interference).

2.4 Working principle

Principle of nebulization:


The mesh nebulizer works by electrically activating a piezoelectric ceramic actuator (PZT), which transmits the generated vibration to the adjacent carrier plate and the polymer mesh with numerous openings. The vibration actively forces the liquid medication out by physically breaking the surface tension of the solution through the mesh openings, resulting in final nebulization. After the patient has put on the mask covering the nose and mouth, the nebulized medication is inhaled into the body through inhalation.

Principle of treatment:

The respiratory system is an open system. After inhalation, the atomized drug can be absorbed directly through the mucous membranes in the patient's oral cavity, throat, trachea, bronchi and alveoli, etc. to achieve the purpose of treatment.

3. Device description

The nebulizer consists of the main unit, the medication cup and the mask.

1. Medication container
2. Nebulizer sheet
3. Mouthpiece
4. Electrodes
5. Main unit
6. On/off button

Accessories:

Adult mask

Children's mask

Mouthpiece

Charging cable

4. Application

4.1 Before initial operation


1. Remove the packaging

Attention: Please clean and disinfect the appliance before using it for the first time.

2. Assembling the inhaler

(1) Attach the medication container to the main unit by sliding it from the front towards the main unit.


Attention: If you attach the medication container to the main device, make sure that you attach it correctly. The medication container is correctly seated when you hear a locking sound. Otherwise, the electrode line may fail. The device will then not be able to atomize normally.

1. Press

(2) Remove the medication container from the main unit. Press and hold the "PUSH" button on the main unit and push the medication cup away from the main unit.

Attention: To prevent damage to the appliance, always press the "PUSH" button before removing the medication container.

1. Push away


2. "PUSH"- Button

4.2 Use

Preparation before use

1. Remove the medication container and clean and disinfect it before use.
2. Adding medication: Open the container lid and add the medication to the medication container as shown below:


1. Container lid
2. Medication insertion


Attention:

- (1) Before using any medicine or medication, please consult your doctor to ensure that you are using the product correctly.
- (2) Do not use medication with a high concentration, high viscosity, oily medication or volatile liquid medication, as this can lead to abnormal atomization.
- (3) It is recommended not to exceed the capacity of the medication container. If the medication container is filled with medication, cover the container lid to prevent leakage.

3. Close the container lid.

4. Attach the medication container to the main device.
5. Attach the mask as shown.

**Use:**

- 1. Switch the device on:** Press the On/Off button. After 1 second, the power indicator lights up green and the device starts to atomize.
- 2. Change mode:** The device has two operating modes: level I and level II. Level I has a high nebulization rate, level II has a low nebulization rate. The operating mode can be changed as required by briefly pressing the on/off button.
- 3. Inhalation:** Hold the device in your hand, place the mask on your face and slowly inhale the medication mist.

Attention:

- If the medication container is not filled with medication, the device switches off automatically. Before switching off, the indicator light flickers orange.
- After pressing the on/off button, the device runs through a short start-up phase (approx. 2 seconds) and then atomizes.
- When the medication runs out, it is recommended to tilt the device slightly towards the user (the side with the button is closer to the user) so that the remaining solution can come into contact with the nebulizer sheet for nebulization.

- Due to the different properties of the medical solutions, the device may not switch off automatically when the medication has been used up, so the user must switch the device off manually to protect the nebulizer sheet from damage.
- Please hold the device steady during use and do not shake it vigorously.
- There is an air hole on the lid of the medicine cup, do not cover it to ensure normal atomization.
- The duration of each inhalation should not exceed 20 minutes. If you experience any discomfort during use, please discontinue treatment and consult a doctor.

4. Switch off the device: If you need to switch off the appliance during use, please press the On/Off button. After 1 second, the indicator light goes out and the appliance switches off.

Attention:

- To prevent the device from being operated accidentally (e.g. when carried in a pocket), the product is designed to be switched on and off after a long press of 1 second.
- After the end of use, it is normal for some medication to remain in the medication container after the automatic switch-off.

4.3 Charging the lithium battery

Operating time of the battery:

- (1) The indicator lights up orange when the battery is low, prompting the user to charge the device.
- (2) In normal use, the battery can work continuously for 1 hour after a full charge.

Charging:

- (1) Plug the charging cable into the port on the device.
- (2) Plug the USB connection of the charging cable into an adapter and plug it into a socket.

During the charging process, the indicator lights up blue and goes out automatically when charging is complete.

Attention!

- **Specification of the power adapter: Input: AC100-240V 50-60Hz, output: DC5V, 1.0A.**
- **Please unplug the mains adapter from the socket after use and avoid connecting the device to the power supply for long periods of time.**
- **If the device is not used for a long time, it should be recharged every six months, which can significantly extend the life of the battery.**
- **The battery cannot be replaced by the user, please contact a local service center or our company if necessary.**
- **If you choose another adapter, it should meet the requirements of IEC60601-1, and its input should be AC100-240V 50-60Hz, and the output should be DC5V, 1.0A.**

5. Maintenance, transportation and storage**5.1 Cleaning and disinfection**

Clean and disinfect the device after each use. If you do not clean the device, atomization will be impaired by drying out and coagulation of medication.

1. **Remove the medication container, accessories and batteries from the main device.**
2. **Open the container lid and dispose of the remaining medication.**
3. **Pour a 75% ethanol solution into the medication container, close the container lid and leave to act for at least 10 minutes. For better disinfection, you can shake gently.**
4. **Immerse the accessories to be disinfected in a container with ethanol solution and close the container. Use a 75% ethanol solution, which should be left to act for 10 minutes or longer.**
5. **Dispose of the disinfectant in the medication container and remove the accessories from the disinfectant. Clean the medication container and accessories several times with fresh water.**
6. **Fill the medication container with clean water, attach it to the main**

unit and let the unit run for 10 minutes to clean the nebulizer sheet.

7. After cleaning, absorb the water with an unused, clean cloth.
8. Wipe the surface of the main unit with 75% medical alcohol and allow it to air dry, or use a clean soft cloth.
9. After completing all the above steps, store the main unit, medication container and accessories in a dry, clean place.

Attention:

- Do not put the medication cup and accessories in boiling water for disinfection, as it could become deformed. Do not put it in a microwave to dry.
- The parts disinfected with disinfectant must be completely cleaned. Otherwise, the residual disinfectant may cause worsening of symptoms.

5.2 Replacing the medication container

The nebulizer sheet is a type of consumable. Generally, the nebulizer life is about six months (at 20 minutes per use, three times a day). Life depends on the application, medication and cleaning. If there is no or little nebulization during use, replace the medication container in a timely manner (please contact your local distributor for assistance).

5.3 Transportation and storage

Transportation and storage environment:

Temperature: -40°C ~ +55°C

Relative humidity: 5 % ~ 96 %

Air pressure: 500 hPa ~ 1060 hPa

Requirements for transportation and storage

- No corrosive gases and a well-ventilated room
- Keep the device out of the reach of children
- Do not store the device where it will be exposed to direct sunlight, high temperatures, humidity, dust or water
- Avoid subjecting the device to tilt, vibration or shock.
- Transportation should be carried out in general means of transportation or in accordance with contractual requirements.

Avoid violent impact, vibration, rain and snow splashes during transportation.

5.4 Environmentally friendly disposal and recycling

The lifespan of the product is 3 years. If the device exceeds the period of use, it must be disposed of. Please contact the manufacturer or distributor for more information.

- 1) The expired nebulizer can be returned to the manufacturer or distributor for proper recycling.
- 2) Used parts can be returned to the manufacturer or distributor for disposal or in accordance with applicable laws and regulations.

6. Troubleshooting

Error	Possible reason	Solution
The device does not start.	Battery is low.	Charge the device.
No or only slight atomization during operation	The medication container is not inserted correctly.	Check the medication container and reinsert it.
	No medication in the container	Add the medication to the container, observing the maximum capacity.
	Inappropriate medication	Ask your doctor whether the medication is suitable for the device.
	The nebulizer sheet is dirty	Clean the medication container.
Water build-up at the nozzle of the atomizer.	Due to temperature differences, the temperature of the surface of the medication cup is relatively low. The medication mist comes into contact with the nozzle and then condenses into water droplets.	Remove the medication container and pour out the water.

After switching on, the power indicator will light up for about 1 second and then go out immediately.	The medication container is not inserted correctly.	Insert the medication container correctly.
	No medication in the medication container	Put medication in the container after consulting your doctor
After switching on, the power indicator lights up 1x, then the device immediately goes off or cannot operate normally.	The battery is dead.	Charge the device.
The nebulizer does not switch off automatically when the medication has been used up.	Medications can cause bubbles in the medication container.	Press the "ON/OFF" button to turn off the device and remove the bubbles.
	Medications can be found on the atomizer sheet.	Press the "ON/OFF" button to turn off the device and clean the medication container.
	The electrodes on the container could be dirty.	Press the "ON/OFF" button to turn off the device and clean the electrodes.
The operating time is too short after the device has been charged.	The battery is not fully charged.	Please charge the device.
	The battery is damaged.	Please contact your local customer service center.
If the device is still not functioning normally after performing all of the above measures, please contact customer service.		

7. Explanation of symbols

Symbol	Description	Symbol	Description
IP22	Degree of water and dust resistance: IP22		Class II device
	Type BF applied part		Keep dry

	Fragile, handle with care		Top
	Humidity for transport: 5% ~ 96%		Temperature range for transportation: -40 °C ~ +55 °C
	Barometric pressure for transport: 500 hPa ~ 1060 hPa		Follow the operating instructions
	Indoor use		On/off button, switch between nebulization rates
	Do not reuse		Follow the instructions for use
	Sterilized using ethylene oxide		Sterile
	Catalog number		Do not use if the packaging is damaged
	Shelf life		Batch number
	Manufacturer		Serial number
	Date of manufacture		Product is compliant with all relevant UK regulations.
	Medical device		Authorized representative in the European Community
	The stacking layer limit is N, depending on the actual conditions.		Disposal according to the Waste Electrical and Electronic Equipment Directive WEEE (2002/96/EC)
	This item complies with the Medical Devices Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community.		

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
Appendix I

Curve chart of equivalent volume particle diameter distribution:

The median particle diameter ($D_{0.50}$) is $1 \sim 4 \mu\text{m}$. Error shall be within $\pm 25\%$

Appendix II Electromagnetic Compatibility (EMC)

Warning

- Don't near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- Use of the NE-M01L adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, NE-M01L and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of NE-M01L could result in increased electromagnetic emissions or decreased electromagnetic immunity of NE-M01L and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the NE-M01L, including cables specified by the manufacturer. Otherwise, degradation of the performance of NE-M01L could result.

Note:

IN-100 (NE-M01L) needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.

The basic performance: Nebulization rate: $\geq 0.25 \text{ mL/min}$.

When the device is disturbed, the data measured may fluctuate, please measure repeatedly or in another environment to ensure its accuracy.

Other devices may affect this device even though they meet the requirements of CISPR.

The following cable types must be used to ensure that they comply with interference radiation and immunity standards:

Name	Cable length(m)
Power cord	1.0

Table 1

Guidance and Declaration - Electromagnetic Emissions	
Emissions test	Compliance
Radiated RF EMISSIONS CISPR 11	Group 1
Radiated RF EMISSIONS CISPR 11	Class A

Harmonic distortion IEC 61000-3-2	Class A
Voltage fluctuations and flicker IEC 61000-3-3	Complies

Table 2

Guidance and Declaration - Electromagnetic Immunity		
Immunity Test	IEC 60601 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines
Surge IEC 61000-4-5	± 0,5 kV, ± 1 kV, line(s) to line(s)	± 0,5 kV, ± 1 kV, line(s) to line(s)
Voltage dips and Voltage interruptions IEC 61000-4-11	0 % UT; 0,5 .cycle .At0°,45°,90°,135°,180°,225°,270°and315°. 0 % UT; 1 cycle and 70 % UT ; 25/30 cycles ;Single phase:at 0°. 0 % UT ; 250/300 cycle	0 % UT; 0,5 .cycle .At0°,45°,90°,135°,180°,225°,270°and315°. 0 % UT; 1 cycle and 70 % UT ; 25/30 cycles ;Single phase:at 0°. 0 % UT ; 250/300 cycl
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V 0,15MHz - 80 MHz 6 V in ISM bands between 0,15MHz to 80 MHz 80%AM at 1kHz	3 V 0,15MHz - 80 MHz 6 V in ISM bands between 0,15MHz to 80 MHz 80%AM at 1kHz
Radiated RF IEC61000-4-3	10V/m 80 MHz-2,7GHz 80%AM at 1kHz	10V/m 80 MHz-2,7GHz 80%AM at 1kHz
NOTE UT is the a.c.mains voltage prior to application of the test level		

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity						
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	IEC60601-1-2 Test level (V/m)	Compliance level(V/m)
	385	380 - 390	TETRA 400	Pulse modulation b) 18 Hz	27	27
	450	430-470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	28	28
	710 745 780	704 - 787	LTE Band 13, 17	Pulse modulation b) 217 Hz	9	9
	810 870	800 - 960	GSM 800/900, TETRA 800,	Pulse	28	28

	930		iDEN 820, CDMA 850, LTE Band 5	modulation b) 18 Hz		
	1720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	28	28
	1845					
	1970					
	2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	28	28
	5240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	9	9
	5500					
	5785					

Novidion GmbH

ADD: Fuggerstr. 30, 51149 Köln, Germany

Tel.: 02203 – 9885 200, Fax: 02203 – 9885 206

info@pulox.de, www.pulox.de

WEEE-Reg.-Nr.: DE24355330

CONTEC MEDICAL SYSTEMS CO, LTD

Add: No.112 Qinhuang West Street, Economic & Technical
Development Zone, Quinhuangdao, Hebei Province
066004, P.R. China

Prolinx GmbH

ADD: Brehmstr. 56, 40239 Duesseldorf, Germany

Tel: 0049 211 3105 4698

E-mail: med@eulinx.eu

Tinovamed GmbH

Add: Zielmattenring 6, 4563 Gerlafingen, Switzerland

Tel: +41 32 566 71 72

E-Mail: info@tinovamed.ch

Web: www.tinovamed.ch

MedPath Limited

Add: 27 Old Gloucester Street, WC1N 3AX London, United Kingdom

Tel: +44 776 280 3783

E-mail: UK@medpath.pro

Web: www.medpath.pro

Bei Fragen wenden Sie sich bitte an:/ If you have any questions, please contact:/ Si tiene alguna pregunta, póngase en contacto con nosotros:/ Si vous avez une question s'il vous plaît contacter:/ Per qualsiasi domanda, si prega di contattare / Als u vragen heeft, kunt u contact opnemen met / W razie jakichkolwiek pytań prosimy o kontakt / Om du har några frågor, vänligen kontakta:

Novidion GmbH, Fuggerstr. 30, 51149 Köln, Germany

Tel.: +49 2203 – 9885 200, Fax: +49 2203 – 9885 206

info@pulox.de, www.pulox.de, WEEE: DE24355330

Sollten Sie beim Lesen dieser Anleitung Probleme haben, so können Sie diese auch auf unserer Website downloaden:/ In case you have any problems reading this manual, you can download it from our website:/ En caso de que tenga algún problema para leer este manual, puede descargarlo de nuestro sitio web:/ Si vous rencontrez des problèmes pour lire ce manuel, vous pouvez le télécharger depuis notre site Web:/ In caso di problemi nella lettura di questo manuale, è possibile scaricarlo dal nostro sito Web / Als u problemen hebt met het lezen van deze instructies, kunt u ze ook downloaden van onze website / Jeśli masz problemy z odczytaniem tej instrukcji, możesz ją również pobrać z naszej strony internetowej / Om du har problem med att läsa dessa instruktioner kan du också ladda ner dem från vår webbplats:

<https://www.pulox.de/Anleitungen-Datenblaetter-Retouren>

Date of the last modification: 2024/12/11 / Version: 1.5