Zuhause TEST

DE Eisenmangel

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EN Iron Deficiency

INTRODUCTION

The **ZuhauseTEST Iron Deficiency** can be used as a tool for the detection of iron deficiency. Anemia due to iron depletion is widely held in children and women of all ages but mainly by women who still have their periods (at least 20% suffer from iron deficiency). Main signs are paleness, feeling tired, headaches, faster heartbeat, or shortness of breath during exercise. They may appear gradually and could go unnoticed. Iron deficiency occurs when blood does not contain enough red blood cells and thus low levels of hemoglobin, which is the major protein involved in oxygen transport in whole body. An important component of hemoglobin is iron.

Depletion of iron, which can happen during pregnancy, growth, in case of insufficient iron intake, inadequate absorption or blood loss (period, abnormal bleedings, ulcers, etc.) has tremendous effects on health. Low ferritin may also indicate hypothyroidism, vitamin C deficiency or celiac disease. Low ferritin levels are seen in some patients with restless legs syndrome, not necessarily related to anemia, but perhaps due to low iron stores just prior to anemia.^{1,2}

The **ZuhauseTEST Iron Deficiency** is a rapid chromatographic immunoassay for the qualitative detection of ferritin in human fingertick blood at a cut-off concentration of 30 ng/mL. The membrane is pre-coated with anti-ferritin polyclonal antibody on the test line region. The gold layer is pre-coated with anti-ferritin monoclonal antibody and Rabbit IgG. During testing, the specimen reacts with the particle coated with anti-ferritin monoclonal antibody. The mixture migrates upward on the membrane chromatographically by capillary action to react with anti-ferritin polyclonal antibody on the membrane and generate a colored line. The line in test line region (T) appears, if the ferritin level exceeds the cut-off level of 30 ng/mL. If the ferritin concentration is less than 30 ng/mL, the test line does not appear. To serve as a procedural control, a colored line will always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

TEST CONTENTS

- 1 Test cassette [Ferritin] with 1 pipette in a sealed pouch
- 1 Tube with sample dilution buffer
- 2 Automatic lancing devices

PROMISEMED HANGZHOU MEDITECH CO., LTD.
No. 1388 Gangping Street, Cangqian Community, Yuhang District, Hangzhou City, 311212 Zhejiang, China

1 Alcohol swab
Jiangsu Sunclan Medical Co., Ltd.
No. 11 Fenghuang south Road, Hutang Town, Wujin District 213162 Changzhou City, Jiangsu Province, P.R. China.

oder
LIGHTS MEDICAL MANUFACTURE CO., LTD.
No. 19 Quanda Road, Wujing Development Area, Tianjin, China

oder
Baoying Country Fukang Medical Appliance Co., Ltd.
Guangyang road, Huangcheng town industrial area, Baoying County, Yangzhou, Jiangsu, 225800, China

Symbolerklärungen | Explanation of symbols | Explication des symboles | Indice dei simboli

 Gebrauchsanweisung beachten Follow instructions Suivre les instructions Consultare le istruzioni per l'uso	 <i>In-vitro</i> -Diagnostikum (Anwendung außerhalb des Körpers) <i>In vitro</i> diagnostic medical device (for external use) Diagnostic <i>in vitro</i> (usage externe) Dispositivo diagnostico <i>in vitro</i> (uso al di fuori del corpo umano)	 Verwendbar bis (siehe Aufdruck Packung) Best before (See imprint on package) À utiliser avant (voir l'impression sur l'emballage) Usare entro (vedere la data riportata sulla confezione)	 Bei Beschädigung nicht verwenden Do not use if package is damaged Ne pas utiliser si l'emballage est endommagé Non usare se la confezione è danneggiata	 Bevollmächtigter in der Europäischen Gemeinschaft Authorized representative in the European Community Représentant autorisé dans la Communauté européenne Rappresentante autorizzato nell'UE	 Bei 2–30 °C trocken lagern. Nicht einfrieren. Store at 2 °C–30 °C. Do not Freeze. Conserver dans un endroit sec entre 2 et 30 °C. Ne pas congeler. Conservare da 2 a 30°C. Non congelare.	 Schweizer Bevollmächtigter Authorised representative for Switzerland Rappresentante autorizzato in Svizzera Représentant autorisé en Suisse
 Hersteller Manufacturer Fabricant Produttore	 Chargenbezeichnung (siehe Aufdruck Packung) Batchnumber (See imprint on package) Numéro de lot (voir l'impression sur l'emballage) Numero di lotto (vedere la data riportata sulla confezione)	 Inhalt ausreichend für 1 Prüfung Content sufficient for 1 test Contenu suffisant pour 1 test Contenuto sufficiente per 1 test	 Bestellnummer Catalogue number Numéro de commande Catalogo #	 Importeur Importer Importateur Importatore	 Nicht wiederverwenden Do not reuse Ne pas réutiliser Non riutilizzare	

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PROMISEMED HANGZHOU MEDITECH CO., LTD.
No. 1388 Gangping Street, Cangqian Community, Yuhang District, Hangzhou City, 311212 Zhejiang, China

1 Compressé d'alcool
Jiangsu Sunclan Medical Co., Ltd.
No. 11 Fenghuang south Road, Hutang Town, Wujin District 213162 Changzhou City, Jiangsu Province, P.R. China.

ou
LIGHTS MEDICAL MANUFACTURE CO., LTD.
No. 19 Quanda Road, Wujing Development Area, Tianjin, China

ou
Baoying Country Fukang Medical Appliance Co., Ltd.
Guangyang road, Huangcheng town industrial area, Baoying County, Yangzhou, Jiangsu, 225800, China

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BIBLIOGRAPHIE

- Kryger MH, Otake K, Foerster J (March 2002). „Low body stores of iron and restless legs syndrome: a correctable cause of insomnia in adolescents and teenagers".Sleep-Med.3(2): 127-32.
- Mizuno S, Mihara T, Miyaoka T, Inagaki T, Horiguchi J (14 March 2005). „CSF iron, ferritin and transferrin levels in restless legs syndrome".J Sleep Res1: 43-75.

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oppure

Baoying Country Fukang Medical Appliance Co., Ltd.
Guangyang road, Huangcheng town industrial area, Baoying County, Yangzhou, Jiangsu, 225800, China

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
EN Wash your hands with soap and rinse them with clear, warm water. Bring the foil pouch containing the test to room temperature before opening it. Open the pouch, take out the test cassette and the pipette, and place them on a clean and flat surface.

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
EN Carefully twist off the cap of the lancing device and dispose of it.

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
EN Use the provided alcohol pad to clean the fingertip of the middle or ring finger as the puncture site. Allow to air dry.

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
EN Press the lancing device against the finger on the side where the cap was removed. After use, the tip will automatically and safely retract.

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
EN Hold your hand slightly downward and massage it, without touching the puncture site, towards the fingertip to obtain a drop of blood.

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
EN Without squeezing the capillary dropper bulb, put it in contact with the blood. The blood migrates into the capillary dropper through the capillarity to the line indicated on the capillary dropper. You may massage again your finger to obtain more blood if the line is not reached. As far as possible, avoid of air bubbles.

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
EN Dispense the **entire blood** sample into the sample well (S) of the cassette by squeezing the bulb of the pipette.

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
EN Wait for the blood to be totally dispensed in the well. Unscrew the cap of the buffer bottle and add **1 drop of buffer** into the Buffer well (S) of the cassette.

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EN Read the results **exactly after 5 minutes**. The result must not be interpreted after 10 minutes.

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EN The test result indicates a **sufficient iron level** if a **colored line** appears in the **control area (C)** and a **colored line** is visible in the **test area (T)**. This result means that the ferritin concentration in the blood is normal and there is no potential iron deficiency. If symptoms indicating a deficiency persist, we recommend consulting a doctor.

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EN The test result indicates an **insufficient iron level** if a **colored line** appears in the **control area (C)** and **no colored line** is visible in the **test area (T)**. This result means that the ferritin concentration in blood is too low. You should consult physician because it may be an iron deficiency.

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EN If you do **not** see a **control line (C)** or **only a test line (T)**, the test did not run correctly and is invalid. Check whether you have followed all steps of the instructions for use precisely. Perform a new test.

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