Opraclean®

Mullbinde zur Tamponade mit Iodoform
Compresse pour mèche avec iodoforme
Ribbon gauze with iodoform for packing
Venda de gasa para el taponamiento con yodoformo
Ligadura de gaze para tamponamento com iodofórmio
Benda di garza per il tamponamento con iodoformio
Gaasverband voor tamponnade, met jodoform
Gazebind til tamponade med iodoform
Gasbinda för tamponad med jodoform
Mulové obinadlo k tamponádě s jodoformem
Mulové ovínadlo na tamponádu s jodoformom
Jodoformos mullpólya tamponáláshoz

Schlauchgaze mit Iodoform
Tube de gaze avec iodoforme
Tubular gauze with iodoform
Gasa tubular con yodoformo
Gaze tubular com iodofórmio
Garza tubolare con iodoformio
Kousgaas met jodoform
Tubegaze med iodoform
Tubgasbinda med jodoform
Tubulární gáza s jodoformem
Tubusová gáza s jodoformom
Jodoformos csőgéz

Verbandmull mit Iodoform
Gaze avec iodoforme
Absorbent cotton gauze with iodoform
Gasa hidrófila con yodoformo
Gaze com iodofórmio
Garza per bendaggi con iodoformio
Verbandgaas met jodoform
Forbindsgaze med iodoform
Gasväv med jodoform
Obvazový mul s jodoformem
Obväzový mul s jodoformom
Jodoformos mullkötszer

25 °C
CE 0123
Do not use if package is damaged or opened
Do not freeze or refrigerate
22 PAP
Recyclable Packaging
4045353/2019-04

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en Instructions for Use

Composition:
100 cm² saturated ribbon gauze for packing contains:
Iodoform: 23 mg
100 cm² saturated tubular gauze contains:
Iodoform: 28 mg
100 cm² saturated absorbent cotton gauze contains:
Iodoform: 57 mg

Other components: cotton, liquid paraffin

Intended purpose:
The Opraclean® dressing materials (ribbon gauze with iodoform for packing, tubular gauze with iodoform, absorbent cotton gauze with iodoform) are used for antiseptic packing and drainage after surgical removal of infected tissue.

Indications:
The Opraclean® dressing materials (ribbon gauze with iodoform for packing, tubular gauze with iodoform, absorbent cotton gauze with iodoform) are applied after surgical removal of infected tissue in:
- abscesses
- fistula-associated abscesses
- ulcers
- wound pockets and body cavities

Contraindications:
- hypersensitivity to iodine
- hypersensitivity to any of the product's other components
- hyperthyroidism
- pregnancy
- careful determination of indication during breastfeeding, in infants, newborns and young children, as well as in thyroid function disorders, in particular in the case of repeated and large-scale applications

Side effects:
Allergic reactions are possible (see Contraindications).
In isolated cases involving large wounds, confusion, hallucinations, loss of consciousness, vomiting, fever or accelerated heart rate with cardiac arrhythmias may occur.

Interactions:
Be aware of possible interactions between Opraclean® Dressing Materials and enzymatic wound treatment agents.

Dosage and method of application:
Unless otherwise prescribed, Opraclean® dressing materials with iodoform should be applied, changed or removed by a doctor only. The doctor will determine the amount to be used. The total amount of iodoform used should not exceed 2 g per dressing change and day. To remove, use scissors to open the primary packaging (= aluminium composite pouch) at the marked spot. This product is intended for single use only. The Opraclean® dressing materials can be combined with these products Vliwaktiv® absorbent activated charcoal dressing, Vliwasorb® superabsorbant wound dressing and Vliwazell® highly absorbent all-purpose dressing.

Duration of application:
The doctor will determine how long the Opraclean® Dressing Materials will remain in place. Generally, the individual dressing should not be left in the wound area for more than 24 hours. Opraclean® Dressing Materials are changed based on the indication and the wound status. When used continuously, the wearing time must not exceed 30 days.

Precautions:
The product must not be used if the packaging has already been opened or is damaged.
When used for large areas in particular, do not apply an occlusion dressing.
When large amounts are used (e.g. in abdominal cavities), the patient might exhibit anaesthetic-like symptoms due to the affinity of iodoform with chloroform. In the event that this occurs, remove the Opraclean® Dressing Materials immediately.
If allergic reactions and the above side effects occur, the dressing material must be removed immediately. A doctor must be consulted if any side effects occur.
The aforementioned side effects usually resolve completely within 3–8 days of discontinuing treatment and they should be treated symptomatically by the doctor where necessary.
Keep the products out of the sight and reach of children! Once the aluminium pouch is opened, the Opraclean® dressing materials must be used immediately.
Opraclean® dressing materials are not sterile and must not be sterilised.

Note:
As a result of recrystallisation of the iodoform, in isolated cases small shiny crystals may form over the course of the storage period (after about 1 to 2 years).

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