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EN **Instructions for use**

Product description
Grassoliind Neutral dressings are non-medicated, sterile, paraffin gauze dressings for wound treatment.

Composition
Lens weave (tulle) cotton gauze fabric impregnated with an anhydrous ointment consisting of:
- White vaseline
- Diglycerol ester of mono- and dicarboxylic fatty acids
- Microcrystalline hydrocarbon wax

Properties and mode of action
The tulle's open weave design allows exudate to pass into a layer of absorbent dressing placed on top of it, so helping to avoid maceration. The ointment keeps the edges of the wound supple, promotes the formation of granulation tissue and epithelium and prevents the dressing from adhering to the wound. Dressings can be changed painlessly without irritating granulation tissue and new epithelium.

Intended purpose
Grassoliind Neutral are sterile, non-medicated single-use ointment dressings for wound treatment and for keeping wound surfaces and edges supple, especially exudating surface wounds (such as extensive abrasions or lacerations as well as burns or scalds, in plastic and cosmetic surgery, in nail removal, craniotomies, etc.) for covering donor sites in skin grafting and for fixing the split-thickness skin graft; ideal also in dermatology as well as for patients with sensitive skin and sensitivity to medications.

The neutral properties of the non-medicated ointment base make Grassoliind Neutral particularly useful in dermatology and for patients with sensitive skin or those sensitive to certain medicinal substances.

Grassoliind Neutral ointment dressings are intended as primary wound covering and are used in combination with a secondary absorbent wound dressing. Grassoliind Neutral is intended for the treatment of human skin by professional users.

Mode of application
Remove the paraffin gauze dressing from the sterile peel pack keeping both product savers in place. If needed, cut the dressing to the size of the wound with sterile scissors. After removing one of the product savers, place this side of the dressing on the wound and remove the other product saver. It is recommended that a sterile, absorbent dressing be placed over the paraffin gauze dressing to absorb wound exudate. Unless otherwise instructed by the doctor, replace Grassoliind Neutral every time the dressing is changed.

The wear time of the dressing depends on the condition of the respective wound and the level of exudate.

Contraindications
Do not use Grassoliind Neutral on patients who may be allergic to any of its ingredients.

Product disposal
In order to minimise the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Special precautions
In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary on these populations groups, this dressing should be used with caution following a clinician's recommendation.

Reusing a single-use medical device is dangerous. Reprocessing device, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.

Before opening the product packaging, inspect it carefully for signs of damage. Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.

Incident reporting
For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices): if during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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Paul Hartmann AG