PORTUGUÊS

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en Intended for FLUSHING ONLY of in-situ PIVCs, PICCs, CVCs, and implanted venous access ports

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0050

Becton, Dickinson and Company, 1 Becton Drive, Franklin Lakes, NJ 07417 USA. Made in Spain. bd.com/symbols-glossary
BECTON Becton Dickinson Ireland Limited, Donore Road, Drogheda Co. Louth, A92 YT52, Ireland
Australia - Imported by: Becton Dickinson Pty Ltd, 66 Waterloo Road North Ryde NSW 2113 Australia
New Zealand - Imported by: Becton Dickinson Limited
148 George Street, Level 10, Wellington 6011 New Zealand
BD, the BD Logo, Luer-Lok and PosiFlush are trademarks of Becton, Dickinson and its affiliates. © 2024 BD. All rights reserved.
SSCP is available on https://eifu.bd.com/#/public/Multilanguage

10000223600/07 500078560 Rev. 2024-08

ENGLISH

INSTRUCTIONS FOR USE

Follow instructions for use, manufacturer guidelines, and institutional regulations.

DESCRIPTION
BD PosiFlush™ SP Syringe is a ready to use medical device (according to regulation (EU) 2017/745 of the European Parliament and of the Council). It is a polypropylene syringe containing sterile and non-pyrogenic isotonic 0.9% sodium chloride solution. The syringe content is guaranteed to be sterile, non-toxic and non-pyrogenic.

INTENDED USE/INDICATIONS FOR USE
BD PosiFlush™ SP Syringes are intended to be used FOR FLUSHING intravenous (IV), peripheral venous catheters (PVC), peripherally inserted central catheters (PICC), central venous catheters (CVC), and implanted venous access ports.

INTENDED USER
BD PosiFlush™ SP Syringe is to be used by healthcare professionals experienced in IV therapy or vascular access and the use of devices for IV therapy. BD PosiFlush™ SP Syringes are manual devices that may be operated by a large range of people with various human characteristics, including hand size and strength.

CONTRAINDICATIONS
- Do not use in patients suffering from hypernatremia and fluid retention when the administration of sodium or chloride ions is clinically detrimental.
- Do not use BD PosiFlush™ SP Syringe if a patient has a known allergy to any of its components, materials or 0.9% sodium chloride solution, which may lead to an allergic response or hypersensitivity reaction.

PERFORMANCE CHARACTERISTICS
- Interoperability with Luer Lock Devices (ANSI).
- All syringe barrels (3 mL, 5 mL and 10 mL) have the same ID, all syringe barrel and plunger rod lengths are the same. This pressure is equivalent for all sizes.

WARNINGS
- Do not use if unit package or content is damaged.
- Do not use if product has been left at freezing temperature.
- Verify the expiration date on the package and syringe label. Do not use if product has expired.
- Do not use if syringe tip cap or stopper is damaged.
- Do not use if solution is cloudy or colored, contains a precipitate, or has any type of suspended particulate matter.
- Do not re-use. Re-use may lead to infection or other illness/injury.
- Small parts are a potential choking hazard for children and adult patients according to your facility protocol.
- Keep out of reach of children.
- Possible complications: undue adverse reactions associated with flushing may include sepsis, infections (localized/systemic), mucocutaneous blocked response, exposure to bloodborne pathogens, air embolism, particulate embolism, blood clots, phlebitis, leakage that may lead to hazardous drug/fluid exposure, extravasation, transitory taste or odor during flushing. Use of contaminated normal saline product may lead to infection and possibly to death.
- Not using aseptic technique and failure to adhere to flushing guidelines may lead to catheter related bloodstream infection and/or related injury or death, catheter failure, catheter related complications such as occlusion, infiltration, extravasation, erythema, swelling or pain.

PRECAUTIONS
- Store at controlled temperature (15-25°C). Excursions permitted to 5-40°C.
- Do not re-sterilize.
- Check with drug manufacturer/instructions for use to ensure compatibility with 0.9% sodium chloride solution prior to use. If 0.9% sodium chloride solution is not compatible, follow the drug manufacturer instructions for flushing practices, or first flush the vascular access device (VAD) with a compatible solution such as 5% dextrose in water to remove traces of the medication in accordance with manufacturer and institution policies.
- Use aseptic technique in accordance with specific medical conditions.
- Treatment needs, age, and weight that may require restricted solution or fluid intake when deciding to flush with 0.9% sodium chloride injection.
- Blunt fill cannulas should not be used to withdraw medication if it will not exceed fluid intake guidelines.
- There are no known incompatibilities of flushing with BD PosiFlush™ SP Syringe in pregnant and lactating women.
- BD PosiFlush™ SP Syringe is not to be used with ISO luer compliant components for intravenous applications.
- Do not use if damaged, open or wet product.
- EU Only: Users should report any serious incident related to the device to the manufacturer and competent authority.

DIRECTIONS FOR USE
- Use aseptic technique. To ensure safe medication preparation and administration, please consult and practice the "7 rights" of medication administration (right patient, right drug, right dose, right time, right route, right reason and right documentation).
- Current practice recommendations are to flush before and after each medication administration, between incompatible medications, and intervals when catheters are not in use. However always consider device manufacturer instructions and your institution policies.
- Follow disinfection protocols as per institution and manufacturer instructions.

Use aseptic technique throughout the procedure.
1. Open pouch and remove syringe.
2. Check that syringe tip cap is in place. Inspect clarity of 0.9% sodium chloride.
3. Depress plunger with tip cap on to release the stopper seal. (Fig. 1)
4. Remove tip cap. Avoid touching the syringe tip to prevent the tip from becoming contaminated. (Fig. 2)
5. Push syringe plunger to expel the air. (Fig. 3)
6. Connect BD PosiFlush™ SP Syringe to vascular access device, taking care that there is no touch contamination of the connection. (Fig. 4) Ensure secure connection.
7. Push syringe plunger to inject the required volume of saline following manufacturer's guidelines. (Fig. 5) Inject at a rate that is safe for vessel overpressure. In case of plunger resistance, it is recommended that the connection is checked.
8. After use, dispose of in accordance with recognized procedure in your institution. (Fig. 6)

COMPOSITION PER UNIT
Polypropylene syringe with BD Luer-Lok™ tip.
Sterile, non-pyrogenic 0.9% sodium chloride solution, not formulated with Natural Rubber Latex. Saline solution: sodium chloride 9 g/l (NaCl 0.9%), anhydrous water for injection, hydrochloric acid for pH adjustment.

SYRINGE CONTENT

3 mL	3 mL ± 0.15 mL
5 mL	5 mL ± 0.15 mL
10 mL	10 mL ± 0.15 mL

SHELF LIFE
The BD PosiFlush™ SP syringe is a pre-filled, single use 0.9% sodium chloride solution. Shelf life is determined to improve clinical effectiveness, stops and time involved in the manual preparation of saline syringes.
- The BD PosiFlush™ SP syringe reduces risk of touch contamination that may occur during manual preparation of saline flush syringes.

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