700 mm

400 mm

5 mm

10 mm

5 mm

Hydrocoll® Classic

Hydrokolloider Wundverband
Hydrocolloid wound dressing
Pansement hydrocolloïde
Hydrocolloïdverband
Medicazione idrocolloidale
Apósito hidrocoloide
Penso hidrocolóide
Υδροκολλοειδή επιθέματα
Hydrokoloidní krytí na rány
Hydrokoloidné krytie
Opatrunek hydrokoloidowy
Hidrokolloid kötözer
Хидроколоидна превръзка за рани
Hidrokoloidna obloga za rane
Pansament cu hidrocoloid pentru plăgi
Hidrokolloid yara örtüsü
水胶质伤口敷料
ضمادة للجروح غروانية مائية

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EN Instructions for use

Product description
Hydrocoll Classic® is a sterile self-adhesive wound dressing consisting of an adhesive hydrocolloid layer covered with a polyurethane film. The product is a non-active medical device. It does not contain any medicinal substances, tissue or cells of human or animal origin or nanomaterials. It is sterilized using irradiation.

Hydrocoll Classic is a dressing with a tapered edge and a thicker core. Hydrocoll Classic Thin is a flexible dressing that conforms well to body contours. Both presentations provide a moist wound healing environment.

Composition
Hydrocoll Classic consists of three materials: hydrocolloid, polyurethane film and release liner. The hydrocolloid layer is composed of sodium carboxymethyl cellulose (CMC) and hot melt adhesive. It absorbs liquids and enables the dressing to adhere to the wound and the surrounding skin. The polyurethane film is permeable to moisture vapour, waterproof, and impermeable to liquids and bacteria.

Properties and mode of action
The product absorbs wound exudate. After the exudate has been absorbed, the hydrocolloid layer transforms into a gel which provides a moist wound healing environment. Since the gel is produced during the application, the dressing does not adhere to the wound. Due to the thin flexible material, Hydrocoll Classic conforms well to body contours. The surface of Hydrocoll Classic is water-repellent.


Intended purpose & indications
Hydrocoll Classic is intended for the management of non- to moderately exuding wounds including injuries to the skin that have breached epidermis or dermis, such as pressure ulcers, leg ulcers, donor sites, postoperative wounds and skin abrasions.
Hydrocoll Classic is indicated for a wide range of non- to moderately exuding chronic and acute wounds, such as pressure ulcers, leg ulcers, donor sites, postoperative wounds and skin abrasions.

Hydrocoll Classic Thin is intended for management of non- to low exuding wounds including injuries to the skin that have breached epidermis or dermis, such as donor sites, postoperative wounds and skin abrasions.
Hydrocoll Classic Thin is indicated for a wide range of non- to low exuding acute wounds, such as donor sites, postoperative wounds, skin abrasions.

Hydrocoll Classic is for short-term use (single use time <7 days and accumulated use time <28 days). This is a professional use only and single-use device.

Intended user
The product must be used by or under the supervision of a healthcare professional.

Mode of application

Clean the wound area in accordance with standard procedures. Ensure the peri-wound skin is dry. Products must be kept clean and sterile before using. Do not use if package is damaged. Select a size of hydrocolloid dressing which overlaps the wound edges by at least 2 - 3 cm. If needed, cut the dressing to size using sterile scissors.
Remove the dressing from the peel pouch, tear off the release liner and apply the adherent side to the wound. Remove the remaining liner and smooth the dressing down onto skin.
For dressing removal, press down on the skin surrounding the dressing and gently lift one of its edges. Repeat until all edges are loosened from the skin surface.
The dressing can be kept on the wound for several days. Change the dressing before it is fully saturated, shows signs of leakage, or as indicated by clinical practice. The same product cannot be used for more than 7 days, and a continuous use cannot exceed 28 days.
Frequent dressing changes on wounds with damaged or delicate skin surrounding the wound are not recommended.

Hydrocoll Classic Secal: remove the release liner and press the dressing into the anal cleft. The tip of the dressing can point upwards or downwards depending on the position and depth of the ulcer. Press the edges on firmly.
Hydrocoll Classic Concave: remove the largest part of the release liner, place the middle part of the dressing over the wound, remove the rest of the release liner, then overlap the wings to fit.

Intended patient population
Hydrocoll Classic is intended for patients which require treatment targeting the indications mentioned in the intended purpose.

Clinical benefit
Relieve patient pain during treatment.

Contraindications
Hydrocoll Classic is contraindicated in any patients who may be allergic to any of its ingredients.

Storage and transport conditions
Transportation conditions should ensure that the product is protected from external pressure, sunlight and rain.

Product disposal
In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

Special precautions
In the absence of available data supporting the use of this dressing on sensitive population groups such as infants, children, pregnant or nursing women, and in the absence of data to the contrary on these population groups, this dressing should be used with caution following a clinician's recommendation.
Reusing a single-use medical device is dangerous. Reprocessing devices, in order to reuse them, may seriously damage their integrity and their performance. Information available on request.
Before opening the product packaging, inspect it carefully for signs of damage. Do not use the product if there are visible signs of damage, such as cracks, channels, pinholes, or nonuniform, torn or incomplete seals.

Side-effects
Very rarely though, there will be some potential risks associated with the use of Hydrocoll Classic, such as:
- Adverse reactions (intolerance to any of its ingredients, erythema, hypergranulation etc.)
- Wound infections
- Risk of maceration (fluid leakage)

Incident reporting
For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices): if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.
Date of revision of the text: 2024-06-26

* HARTMANN Hydrocoll Classic is registered as Hydrocolloid Dressing with Tapered Edges (Hydrocoll Classic) and Hydrocolloid Dressing Standard (Hydrocoll Classic Thin) under legal Manufacturer Winner Medical Co., Ltd.
Any reference to Hydrocoll Classic contained in the text refers to all presentations – unless otherwise specified.

GB – PAUL HARTMANN Pty. Ltd. · Macquarie Park, NSW 2113
GR – PAUL HARTMANN Ltd. · Heywood Greater Manchester OL10 2TT
ZA – HARTMANN South Africa · Northriding, 2162

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